


# Lung conditions at-a-glance, using ASV

In ASV mode, the Dynamic Lung and ASV Graph, together with RCexp value, provide visual insight to the patient's lung condition.

## Normal lungs




RCexp (s)	Rinsp (cmH2O s/l)	Cstat (ml/cmH2O)
-----------	-------------------	------------------

0.50 – 0.70	10 – 15	45 – 65
-------------	---------	---------



Wide square-shaped safety window

## Obstructive lung diseases



RCexp (s)	Rinsp (cmH2O s/l)	Cstat (ml/cmH2O)
-----------	-------------------	------------------


> 0.70*	16 – 33	50 – 80
---------	---------	---------



Narrow and high safety window

\* A long RCexp (> 0.70 s) indicates increased resistance due to the patient and/or the endo-tracheal tube: COPD, asthma, bronchospasm, endotracheal tube obstruction, and the like.<sup>1</sup>

## Restrictive lung diseases



**RCexp (s)**      **Rinsp (cmH2O s/l)**      **Cstat (ml/cmH2O)**


.....  
 < 0.50\*      10 – 15      < 45

\*A short RCexp (< 0.50 s) indicates decreased compliance due to the lung and/or the chest wall: ARDS, lung fibrosis, atelectasis, kyphoscoliosis, increased abdominal pressure, and the like.<sup>1</sup>



Low and wide safety window  
 (low compliance = stiff lung)

## Mixed conditions



**RCexp (s)**      **Rinsp (cmH2O s/l)**      **Cstat (ml/cmH2O)**

.....  
 0.50 – 0.70      > 16      < 45



Narrow rectangular  
 safety window

1. *Parameters for Simulation of Adult Subjects During Mechanical Ventilation.* Jean-Michel Arnal, Aude Garnero, Mathieu Saoli and Robert L Chatburn  
 Respiratory Care February 2018, 63 (2) 158-168; DOI: <https://doi.org/10.4187/respcare.05775>