

Anesthetic gas delivery systems

Compatibility with Hamilton Medical ventilators

The use of anesthetic gas delivery systems in combination with Hamilton Medical ventilators has been validated successfully as per the following table. Please follow the instructions for use in the operator's manual provided by the manufacturer of the respective product as shown below.

Manufacturer	System	Anesthetic gas removal	HAMILTON-S1	HAMILTON-G5	HAMILTON-C1	HAMILTON-T1	HAMILTON-MR1	HAMILTON-C2	HAMILTON-C3	HAMILTON-C6
Sedana	AnaConDa ^{§1)}	FlurAbsorb 1-day FlurAbsorb 5-day	✓	✓	✓	✓	✓	✓	✓	✓
TIM	Mirus [§]	ORS	✓	✓	✓	✓	✓	✓	✓	✓

Compatible: ✓ Not compatible: X

1) Please note that use of the AnaConDa system may cause a slight increase in airway resistance.

For all proprietary trademarks (®) and third-party trademarks (§) used by Hamilton Medical AG see www.hamilton-medical.com/trademarks.
© 2021 Hamilton Medical AG. All rights reserved.

ELO20210121N.00