

Trolley for the HAMILTON-C2/C3

Transport and stationary solution with a small footprint

The trolley for the HAMILTON-C2 and HAMILTON-C3 ventilators has been developed with a focus on combining a modern design with usability. By virtue of the small footprint, the trolley requires minimal space at the bedside. This makes it ideally suited for intensive and intermediate care units, as well as for intrahospital transport.

Thanks to the design with a low center of gravity, the trolley offers high stability during transport and while stationary. It can also be easily maneuvered over doorway thresholds or other rises.

Features at a glance

- ✓ Small footprint, requiring little space at the bedside
- ✓ Four independent wheels, each with its own brake, for increased maneuverability
- ✓ Brakes can be operated comfortably by foot
- ✓ Trolley column with front and back grooves for mounting accessories
- ✓ Quick-locking mechanism for easy attachment and removal of the ventilator
- ✓ Resistant to cleaning, disinfection, corrosion, and UV light



Swiss
Quality

HAMILTON
MEDICAL

Intelligent Ventilation since 1983

Specific accessories

The following accessories are designed specifically for the new trolley. Previous versions are not compatible.

1 Humidifier holder

The humidifier holder enables the trolley to be fitted with the HAMILTON-H900 or another humidifier and can be adjusted in height (30–60 cm).

2 O2 cylinder holder

The height-adjustable cylinder holder can hold up to two E-type cylinders simultaneously (max. 8 kg for each cylinder) by means of two straps.

Previous compatible accessories

3 Tubing support arm with quick-lock

5 Quick-lock for tubing support arm

4 Water bottle holder (not shown in picture)



The trolley and its accessories can be ordered using the following numbers:

New standard trolley	PN 160170
Humidifier holder	PN 160091
O2 cylinder holder	PN 160090
Tubing support arm with quick-lock	PN 160153
Water bottle holder	PN 160162
Quick-lock for tubing support arm	PN 160154

Manufacturer:

Hamilton Medical AG

Via Crusch 8, 7402 Bonaduz, Switzerland

+41 (0)58 610 10 20

info@hamilton-medical.net

www.hamilton-medical.com